

# Auto-deposition KIT – AC-05-ADK

## Instruction for Use



**GBIEU** (1) Plastic rotative unit, (2) Flat circular magnet ring, (3) Disk, (4) Ring, (5) Clamping ring, (6) Retaining clip



**GBIEU** Insert (2) Flat circular magnet ring into (1) Plastic rotative unit then the (3) Disk and the (4) Ring



**GBIEU** Tighten the (5) Clamping ring in the (1) Plastic rotative unit



**GBIEU** The system is ready for deposition in a beaker on stirring plate



**GBIEU** Place the (6) Retaining clip in the (5) Clamping ring

